

V9.5 Modifications List

Major Features:

- **Extensive User Interface Enhancements** – Many of the major windows and dialog boxes within SpringCharts have been redesigned to increase efficiency and ease of use. Context sensitive popup menus have been added to many items of the program to reduce typing and increase the speed of charting.
- **Enhanced Print Formatting** – V9.5 offers new print formatting options for letters, reports, lists and more. Users can now select font type, font size and margins to customize printed output to meet their needs.
- **Digital Forms** – New forms feature allows users to create electronic forms such as patient intake forms, consent forms and NPP forms. These forms can be completed on-screen, digitally signed and stored within a patient's chart.
- **Digital Signature Support** – Users can now create a digital signature within SpringCharts and attach it to forms, prescriptions, and other documents.
- **Hosted Practice Solution** – Tested configuration enables multiple copies of SpringCharts to be installed and run on the same server for multiple patient database applications.

Additional New Features/Updates:

Features:

- Added the capability to search for drugs by generic name
- Immunizations can now be sorted and printed by type and date
- Out-of-compliance chart evaluations display when a chart is opened and can be processed easily
- Added NPI number to Routing Slips
- Added Lock icon to signed encounters in the Caretree, enabling locked encounters to be identified at a glance
- Addendums can now be added to charted Messages
- Tests have been added to OV Templates
- PopUp Context Menus have been added to text areas (right click to copy/paste/cut, etc.)
- PopUp Context Menus have been added to vitals panels (right click to activate drop down lists of numbers)
- PopUp Context Menus have been added to each category (patient demographics, problem list, insurance, etc.) in the Patient Overview Panel (right click to activate, functions vary by category)
- One encounter type can now be changed to another (for instance, a Vitals Only encounter can be changed to a complete Office Visit)
- A carriage return has been added between Chief Complaint and Present Illness
- PMS import functions no longer require a patient's Social Security Number
- Added a function to delete all routing slips
- Changed 'SOAP' headings (letter only) to Subjective, Objective, Assessment, Plan headings (entire word) on Office Visit GUI and printed Office Visit (i.e. changed 'O' to Objective)
- Midnight shutdown requirement has been removed. Instead, the Practice Management window automatically resets (tracker, appts, todos, etc.) at midnight

- Category Preferences can now be reordered and the dialog is tabbed for easy navigation
- A diagnosis button has been added to the Transactions dialog that enables previous diagnoses for the patient to be added to a transaction
- A confirmation is now received when a report is emailed to a person selected from the address book

User Interface Updates:

- Window and dialog redesigns – Most of the major screens in SpringCharts have been redesigned to provide additional information and/or expedite common functions including:
 - Office Visit
 - Office Visit Template Editor
 - Face Sheet
 - Nurse Note
 - Rx Refill
 - TC Note
 - Patient Chart Overview Panel
- “Referring Dr” and “Referred By” have been changed to “Other Providers” in most areas
- Routing slip text now prints in a smaller font size
- A user preference has been added for right/left orientation for Office Visit window
- The Spell Checker window has been refined

SpringLabs Changes:

- Added 'H' (high), 'L' (low) and 'A' (abnormal) flags to charted lab results to provide at-a-glance indications of out-of-range items
- In lab result reports, if a normal value is empty, no abnormal flag appears
- Changed 'Reference Lab' menu select to 'Reference Lab/Reports'
- Changed timing for standard lab import to a 15 second delay (was 50 seconds)
- Reference/Pending/Completed Lab dialogs now have column-based sort features
- Abnormal Highlights now appear in String Type Lab Items (like Pap smears) in the 'Completed Lab' dialog
- Test Items can now be reordered in the 'New Lab Test' dialog
- PopUp Text has been added to 'Notes' sections of:
 - Completed Lab Test GUI
 - Completed Medical Test GUI
 - Completed Imaging Test GUI

Fixes:

- Fixed transaction message for non-Medical Connect transactions
- Fixed problem of doctors not deleting from user list correctly
- Fixed an issue with patient manager so that patient updates are only broadcast to relevant connected clients when a patient demographic is modified (via SC or PMS)
- Fixed CCR export for Med Tests
- Fixed Chart Evaluations for legacy (pre-version 9) systems.
- Fixed report header for test reports

- Corrected a problem with dialogs not closing with Java 6
- PopUp Text now opens to the correct area when editing
- Changed Routing Slip export data to match the data in printed Routing Slips.
- Duplicate diagnoses are now prevented in Office Visit
- Logout timer now resets on clicks and keyboard actions
- Fixed Imaging and Med Test dialogs – window did not disappear when user clicked 'save'
- Growth Chart Screens have been resized so that the scroll pane shows all of the data